

# Immunocytochemistry of cultured cells

## *Solutions*

### **Paraformaldehyde fixation solution (PFA)**

4% paraformaldehyde

PBS (use TBS instead of PBS hereafter when an alkaline phosphatase is the marker enzyme)

*Heat the solution ( $\leq 70^{\circ}\text{C}$ ) to dissolve paraformaldehyde. Note: Paraformaldehyde is hazardous!*

### **Methanol or Acetone fixation**

Ice-cold methanol

Ice-cold acetone

### **Blocking solution**

1-3% Bovine Serum Albumin (and 10% goat serum)

0.1-0.3% Triton X100 or 0.5% TWEEN 20 (detergent is used when cells are fixed with PFA)

0.05% Sodium Azide ( $\text{NaN}_3$ , Sigma S-2002)

PBS

### **Antibody (Ab) solutions**

Ab-solutions with desired concentrations of primary and secondary Ab are mixed in blocking solution. Add 250  $\mu\text{l}$  Ab-solution per well in a 24 well cassette.

## *Procedure*

If the cells are fixed with PFA use a detergent in the blocking solution from step 3 to 6. If cells are fixed with methanol or acetone the detergent can be omitted. Prepare coverslips with only the secondary Ab to check for unspecific immuno staining.

1. Wash with PBS or  $\text{Ca}^{2+}$  medium 2-3 times. Use  $\text{Ca}^{2+}$  medium if cells are pre-treated.
2. Fix cells with ice-cold PFA, methanol, or acetone fixation solution for 10 min.
3. Incubate with blocking solution for 30-60 min, shaking at room temperature (RT).
4. Incubate with primary Ab-solution 1 h to overnight, shaking at  $4^{\circ}\text{C}$ .
5. Wash with blocking solution 3 times (1 $\rightarrow$ 5 $\rightarrow$ 15 min), shaking at RT.
6. Incubate with secondary Ab-solution under aluminum foil for 1 h, shaking at RT.
7. Rinse with clean PBS 5 times (1 $\rightarrow$ 5 $\rightarrow$ 15 $\rightarrow$ 30 $\rightarrow$ 60 min) under aluminum foil, shaking at RT.
8. Mount the coverslips using Prolong Antifade Kit (Molecular Probes No. P-7481).
9. Leave the slides in dark RT overnight to dry the mounting medium.
10. Examine slides with microscope within a couple of days and then store slides in  $-20^{\circ}\text{C}$ .

Unspecific background fluorescence can be reduced if using filtered blocking solution in step 5 and 6. Filter the solution with a syringe driven filter unit 0.22  $\mu\text{m}$  (MILLEX-GS 33 mm).